

# A Headspace Gas Detection Tube Method to Measure SO<sub>2</sub> in Wine without Disrupting SO<sub>2</sub> Equilibria

Jussara M. Coelho,<sup>1,4</sup> Patricia A. Howe,<sup>2,5\*</sup> and Gavin L. Sacks<sup>3</sup>

**Abstract:** The headspace gas detection tube (HS-GDT) method to measure molecular and free sulfur dioxide (SO<sub>2</sub>) in wine is a simple and inexpensive procedure using commercial industrial safety colorimetric tubes. A syringe is used to sample a wine and to create a closed headspace, which is expelled through the GDT after equilibrium is obtained. The vapor-phase concentration of SO<sub>2</sub> (P<sub>SO<sub>2</sub></sub>) is determined from the manufacturer's printed markings and then related to the molecular SO<sub>2</sub> concentration in the wine based on Henry's law coefficients or, more accurately, on calibration curves. Typical wine ethanol concentrations had a significant effect on the pK<sub>a</sub> of SO<sub>2</sub>, as previously reported, but no effect on Henry's law coefficients. Calibration curves in model wine and aqueous buffer yielded indistinguishable results when appropriate pK<sub>a</sub> values were used. Best results for calibration curves were achieved with 200 mL headspace of model wines (5 to 40 mg/L free SO<sub>2</sub>, equivalent to 0.14 to 1.12 mg/L molecular SO<sub>2</sub>), which yielded satisfactory linearity (r<sup>2</sup> = .99), reproducibility (mean CV = 8% for molecular SO<sub>2</sub> >0.4 mg/L), and detection limits (0.21 mg/L molecular SO<sub>2</sub>). Molecular SO<sub>2</sub> measured by HS-GDT and by aeration-oxidation (A-O) showed a strong correlation between methods for white and blush wines (r<sup>2</sup> = 0.97) and a poor correlation for red wines (r<sup>2</sup> = 0.72). A-O values averaged double the HS-GDT values in red wines. The difference in molecular SO<sub>2</sub> values by A-O and HS-GDT correlated well with the estimated molecular SO<sub>2</sub> decrease due to formation of anthocyanin-bisulfite adduct (r<sup>2</sup> = 0.936), supporting the hypothesis that dissociation of anthocyanin-bisulfite adducts occurs during A-O analysis. By not perturbing this equilibrium, HS-GDT accurately reports free and molecular SO<sub>2</sub> values.

**Key words:** sulfur dioxide, wine analysis, headspace, method validation, SO<sub>2</sub>

Sulfur dioxide (SO<sub>2</sub>) is widely used as a preservative in winemaking due to its antimicrobial and antioxidant properties (Boulton et al. 1996, Ribéreau-Gayon et al. 2007). At wine pH (3 to 4), the majority of SO<sub>2</sub> (>95%) exists in the form of bisulfite (HSO<sub>3</sub><sup>-</sup>), with the remainder existing as neutral "molecular" SO<sub>2</sub> and a negligible portion existing as sulfite (SO<sub>3</sub><sup>2-</sup>) (Zoecklein et al. 1995, Boulton et al. 1996). The sum of these forms (molecular SO<sub>2</sub>, HSO<sub>3</sub><sup>-</sup>, and SO<sub>3</sub><sup>2-</sup>) is referred to as "free SO<sub>2</sub>". Additionally, HSO<sub>3</sub><sup>-</sup> may form covalent adducts with electrophilic compounds in wine, including ketonic acids, sugars, quinones, and anthocyanins (Burroughs and Sparks 1973, Beech et al. 1979). This covalently bound HSO<sub>3</sub><sup>-</sup> is referred to as "bound SO<sub>2</sub>" and has lower antimicrobial and antioxidant activity than free SO<sub>2</sub> forms (Boulton et

al. 1996). "Total SO<sub>2</sub>" is the sum of the free and bound SO<sub>2</sub> and is the form subject to regulatory limitations in most wine-producing countries (Zoecklein et al. 1995).

Of the forms described, the molecular SO<sub>2</sub> and HSO<sub>3</sub><sup>-</sup> species are reportedly the most important for the germicidal and antioxidant properties of SO<sub>2</sub>, respectively (Danilewicz 2011, Divol et al. 2012). For dry table wines, winemakers will typically aim to maintain free SO<sub>2</sub> (predominantly HSO<sub>3</sub><sup>-</sup> for reasons described earlier) at above 10 mg/L (Godden et al. 2001) or at a target of 30 mg/L (Waterhouse and Elias 2010) to mitigate or slow the appearance of oxidative effects such as aldehyde formation and browning. Recommended molecular SO<sub>2</sub> concentrations in wine production depend on the target microorganism, but the most commonly cited values range from 0.5 to 0.8 mg/L molecular SO<sub>2</sub> (Beech et al. 1979, Boulton et al. 1996). Free SO<sub>2</sub> concentrations can decrease during wine production and storage through binding, volatilization, or, more often, oxidation reactions. Thus, winemakers must routinely assess wines for free SO<sub>2</sub> to ensure stability (Zoecklein et al. 1995).

In most production wineries and research laboratories, free SO<sub>2</sub> rather than molecular SO<sub>2</sub> is measured, and the molecular SO<sub>2</sub> is calculated via a modified Henderson-Hasselbalch equation following measurement of pH. Common approaches for measuring free SO<sub>2</sub> fall into two categories: titrimetric and colorimetric methods that rely on direct addition of an oxidizing reagent to the wine, e.g., by iodometry (Ripper method) or addition of *p*-rosaniline (Joslyn 1955, Zoecklein et al. 1995) or alternatively, methods that separate free SO<sub>2</sub> from the wine, followed by quantification by titrimetry, colorimetry, or other means (Zoecklein et al. 1995, Pundir and Rawal 2013). In modestly equipped wineries, the most common variant of

<sup>1</sup>Department of Food Technology, Federal University of Viçosa, Viçosa, MG, Brazil; <sup>2</sup>Department of Food Science, <sup>3</sup>Associate Professor, Department of Food Science, Cornell University, Ithaca, NY 14853; <sup>4</sup>present address, Professor, Departamento de Engenharia de Alimentos, Centro de Ciências Agrárias, Universidade Federal do Espírito Santo, Alto Universitário, s/n, Bairro Guararema, Alegre/ES; and <sup>5</sup>present address, Sr. Director, Quality, Constellation Brands, 7801 St. Helena Highway, Oakville, CA 94562.

\*Corresponding author (pah237@cornell.edu; tel: 707-294-7882)

Acknowledgments: The authors gratefully acknowledge ETS laboratories for providing analysis on the wine samples. The Brazilian CAPES Foundation funded scholarship support for Jussara Coelho (Process number 18868-12-6). Additional funding was from the Peter and Tacie Saltonstall Endowment and the New York Wine and Grape Foundation.

Supplemental data is freely available with the online version of this article at [www.ajevonline.org](http://www.ajevonline.org).

Manuscript submitted Nov 2014, revised Feb 2015, accepted Feb 2015

Copyright © 2015 by the American Society for Enology and Viticulture. All rights reserved.

doi: 10.5344/ajev.2015.14125

this “separation-first” approach is the aeration-oxidation (A-O) method, in which an acidified wine sample is stripped by a stream of air, from which SO<sub>2</sub> is captured by a hydrogen peroxide solution to generate H<sub>2</sub>SO<sub>4</sub>, which can be quantified by titration (Rankine and Pocock 1970). Automated variants of both direct and separation-first methods also exist (Ruiz-Capillas and Jiménez-Colmenero 2009, Pundir and Rawal 2013).

Beyond these methods, the Association of Analytical Communities (AOAC) lists several methods to determine total sulfites in wine (Monier-Williams, flow injection/chromatography-spectrophotometry, and ion exclusion chromatography) and a method for free sulfites for white wines only (flow injection/chromatography-spectrophotometry) (AOAC 2015). These methods give comparable values to A-O, and the A-O method is the one used by agencies such as the United States Alcohol and Tobacco Tax and Trade Bureau due to its relative simplicity (<http://www.ttb.gov/ssd/pdf/tm500.pdf>).

A well-known shortcoming of the direct titrimetric or colorimetric approaches is that other reducing compounds than free SO<sub>2</sub> are also measured, resulting in overestimates of the true value (Buechsenstein and Ough 1978). However, a less-appreciated shortcoming of all standard methods for SO<sub>2</sub> measurement is the dissociation of weak SO<sub>2</sub> adducts following changes to wine equilibrium during the measurement. These perturbations to equilibrium may result from acidification, temperature change, dilution, and/or consumption of free SO<sub>2</sub> by reagents over the measured timecourse, which will overestimate free and molecular SO<sub>2</sub> (Rankine and Pocock 1970, Burroughs 1975, Boulton et al. 1996). This overestimation is more pronounced in red wines, likely due to dissociation of anthocyanin-bisulfite adducts; more modest discrepancies are observed for white wines (Burroughs 1975, Bogren 1996).

This overestimation of free and molecular SO<sub>2</sub> in wine by standard techniques is of potential concern to researchers or winemakers interested in a mechanistic understanding of oxidation or microbial growth in wine (Boulton et al. 1996). However, few papers report using techniques that avoid pH shifts, sample dilution, and/or temperature changes, and thus avoid disturbance of SO<sub>2</sub> equilibria in wine. One such approach is capillary electrophoresis (CE), which can quantify bisulfite with minimal contributions from artifactual dissociation of bound SO<sub>2</sub> forms (Boulton et al. 1996, Collins and Boulton 1996). Results from CE indicate that free SO<sub>2</sub> may be overestimated by up to an order of magnitude in red wines. Similar observations were reported using a colorimetric method that compared absorbance at 520 nm for acetaldehyde- and SO<sub>2</sub>-treated wines to an original wine sample, although this method relied on several assumptions regarding anthocyanin spectral behavior, was developed before an appreciation of copigmentation effects had been described, and would only be appropriate for detecting overestimation due to anthocyanins and not to other weak binders (Burroughs 1975).

Alternatively, the headspace SO<sub>2</sub> concentration of an equilibrated sample can be measured and the partial pressure related to the molecular SO<sub>2</sub> concentration in wine via Henry's Law. Using headspace gas chromatography (HS-GC) coupled to an electrolytic conductivity detector (Davis et al.

1983), good agreement was observed between A-O and HS-GC approaches for a white wine and juice, but A-O yielded 45% higher values for the red wine. Headspace infrared (HS-IR) spectroscopy can also reportedly be used to quantify SO<sub>2</sub> in unadjusted wine samples (Henningsen and Hald 2003). Inductively-coupled plasma optical-emission spectrometry (ICP-OES; Čmelík et al. 2005) and flame molecular absorption spectrometry (Huang et al. 2008) have also been used to measure vapor-phase SO<sub>2</sub> in wine, but wine samples were acidified and gas-purged during analysis, which was expected to yield results comparable to classic methods. None of these approaches to measuring vapor-phase SO<sub>2</sub> have been widely adopted in either wineries or research labs, possibly due to the cost and complexity of the techniques.

An unexplored approach to measuring wine headspace SO<sub>2</sub> concentrations is colorimetric gas detection tubes (GDTs). GDTs were originally developed for industrial safety in the early 1900's (Haag 2001) and consist of a glass tube packed with a color-sensitive reagent, such that the length of tube undergoing color change is proportional to the analyte concentration. Commercial SO<sub>2</sub> gas detection tubes have been used to measure free SO<sub>2</sub> in wines using a modified A-O protocol (Pegram et al. 2013) adapted from a related method for quantifying H<sub>2</sub>S in wines (Park 2008, Ugliano and Henschke 2010). Similar to A-O, this earlier use of GDT to quantify free SO<sub>2</sub> relied on acidification of the sample and sparging prior to SO<sub>2</sub> measurement, and thus likely also experienced dissociation of weak adducts and overestimated free SO<sub>2</sub>.

In this paper, we report on the development of a headspace gas detection tube (HS-GDT) method to determine molecular and free SO<sub>2</sub>, both calculated from the measurement of the headspace gas concentration. This method does not require sample preparation and thus avoids disturbing SO<sub>2</sub> equilibria, and in combination with its low cost and simplicity, should facilitate investigations into the antimicrobial and antioxidant properties of SO<sub>2</sub>.

## Materials and Methods

**Chemicals.** Potassium metabisulfite (97% w/w) and ethanol (95% v/v) were obtained from Acros Organics (Geel, Germany). Lactic acid (85% w/w), potassium bitartrate (99% w/v), hydrogen peroxide (30% w/v), sodium hydroxide (0.01 N), and *o*-phosphoric acid (85% w/w) were obtained from Fisher Scientific (Waltham, MA). A nominally 25% phosphoric acid solution was prepared as a 2.38:1 dilution of 294 mL phosphoric acid (85%) with 700 mL deionized water. Potassium acetate (99% w/w) and potassium sorbate (98% w/w) were obtained from J.T. Baker (Phillipsburg, NJ). Hydrochloric acid (36.5% w/w) was obtained from BDH Merck (Poole Dorset, UK).

**SO<sub>2</sub> working standards.** SO<sub>2</sub> stock solutions at nominal concentrations of 1000 mg/L as SO<sub>2</sub> were prepared weekly by dissolving potassium metabisulfite in a 10% (v/v) solution of methanol in water to avoid SO<sub>2</sub> autooxidation. Working standards were then prepared as needed by adding an appropriate volume of a stock SO<sub>2</sub> solution to saturated potassium bitartrate buffer (pH 3.58). Iodometric titrations were

used to determine the actual concentrations of the stock and working solutions.

**SO<sub>2</sub> gas detection tubes.** Three commercial SO<sub>2</sub> detection tubes were evaluated initially: i) Gastec 5Lb tubes (Gastec Corporation, Fukayanaka, Japan), ii) Sensidyne/Kitigawa 103SE SO<sub>2</sub> tubes (Sensidyne LP, St. Petersburg, FL), and iii) Draeger 0.1/a SO<sub>2</sub> tubes (Drägerwerk AG & Co. KGaA, Lübeck, Germany).

**SO<sub>2</sub> measurements by A-O.** A-O was performed according to a protocol described elsewhere (Iland et al. 2004).

**Protocol for SO<sub>2</sub> measurement by HS-GDT.** The HS-GDT apparatus was assembled as depicted (Figure 1). It consists of a 60 mL Becton Dickinson syringe (polypropylene barrel, polypropylene plunger, latex-free polyisoprene rubber tip treated with Dow 360 medical-grade silicone) with a Luer-lock tip, fitted to a two-way polycarbonate male Luer stopcock. A short piece of silicone tubing connects a GDT to the top of the stopcock. The syringe is fitted with a 3 cm plastic dispensing stop, used to prevent accidental expulsion of liquid during syringe depression, made by cutting a second 60 mL Becton Dickinson syringe barrel to the correct height. If the GDT had been previously used and the “end” marking from previous runs was unclear or had shifted, the “start” point of color transition was marked on the tube with a fine-point permanent marker.

For each analysis, the syringe was used to aspirate 10 mL wine or working standard. The syringe was inverted and then withdrawn further to create 50 mL headspace. The stopcock was closed, the syringe was placed nose-up, the syringe stop was put in place, and the sample was allowed to equilibrate for five min. The stopcock was then connected to a GDT via the short piece of silicone tubing, the stopcock was opened, and the syringe was depressed at a constant rate to the stop, such that 50 mL of headspace was expelled through the GDT over 10 sec. Following gas expulsion, the colorimetric reaction was allowed to stabilize for one min. For some calibration curve experiments and for analyses of wines, this process

was repeated up to four times for one sample, resulting in up to 200 mL headspace through the GDT. For the wine analyses, the limit for the number of repetitions was determined by the number that could be quantified on one tube. The final color transition point was marked with a fine-tip marker and the locations of the stain startpoint and endpoint were measured in millimeters with a ruler.

**Determination of P<sub>SO<sub>2</sub></sub>, molecular SO<sub>2</sub>, and free SO<sub>2</sub> from raw GDT measurements.** Gas detection tubes have printed scales showing the vapor-phase concentration of



**Figure 1** The headspace gas detection tube apparatus prior to (left) and after (right) headspace expulsion through the gas detection tube.

SO<sub>2</sub> (P<sub>SO<sub>2</sub></sub>) in ppm (μL/L) as a non-linear function of stain distance, based on a brand-dependent standard volume (200 mL for the Gastec and Sensidyne/Kitigawa tubes, 100 mL for the Draeger) of headspace sampled. The scale resolution was insufficient for repeated high-precision measurements. To facilitate interpolation, a cubic function, f(x), was determined to relate SO<sub>2</sub> vapor pressure (P<sub>SO<sub>2</sub></sub>) based on the manufacturer’s markings as a function of distance along the GDT in mm. The best-fit equations were calculated using Microsoft Excel (Microsoft, Redmond, WA). Because the spacing of markings varied among GDT lots, cubic interpolation functions were determined for each lot of tubes. For each sample, P<sub>SO<sub>2</sub></sub> was calculated from the location of the startpoint (x<sub>start</sub>) and the endpoint (x<sub>end</sub>) measurements in mm and adjusted based on the actual volume of headspace gas sampled, shown here for the Gastec or Sensidyne tubes (Equation 1).

$$P_{SO_2} = (f(x_{end}) - f(x_{start})) \times \frac{200 \text{ mL}}{\text{headspace sample (mL)}} \quad \text{Eq. 1}$$

Molecular SO<sub>2</sub> (mg/L) was calculated from P<sub>SO<sub>2</sub></sub> using Henry’s Law and the molecular mass of SO<sub>2</sub> (Equation 2). A Henry’s Law coefficient (K<sub>H</sub>) of 0.38 Atm/M at 294 K was determined as part of this work.

$$\text{molecular SO}_2 \text{ (mg/L)} = \frac{P_{SO_2}}{K_H} \times 0.064 \quad \text{Eq. 2}$$

A literature value of 3100 K was used as the temperature dependence constant (Sander 2015) from which to calculate the temperature-corrected Henry’s Law coefficient, K<sub>H(T)</sub>, used in this work, where T is the temperature during the analysis (Equation 3).

$$K_{H(T)} = 0.38 \text{ Atm/M} \times \exp\left(3100 \times \left(\frac{1}{294} - \frac{1}{T}\right)\right) \quad \text{Eq. 3}$$

Free SO<sub>2</sub> was calculated from molecular SO<sub>2</sub>, the sample pH, and the pK<sub>a</sub> of SO<sub>2</sub> (Equation 4). Literature values for pK<sub>a</sub> as a function of ethanol were used in this calculation (Usseglio-Tomasset and Bosia 1984) and were corroborated as part of this work.

$$[\text{free SO}_2] = (1 + 10^{(pH - pK_a)})[\text{molecular SO}_2] \quad \text{Eq. 4}$$

**Evaluation of interferences.** The potential of ethanol, lactic acid, sorbic acid, and acetic acid to act as interferences were evaluated by preparing solutions containing either 14% ethanol, 5 g/L lactic acid, 300 mg/L sorbic acid, or 1.4 g/L acetic acid in a potassium bitartrate buffer. The HS-GDT responses for these solutions were compared against a blank potassium bitartrate buffer. Analyses were done in triplicate using 50 mL headspace samples. For these evaluations, a measurable discoloration of the GDT was taken as evidence of interference.

**Determining optimal equilibration time.** To determine the minimum time required to establish vapor-liquid equilibrium of SO<sub>2</sub>, a 50 mg/L SO<sub>2</sub> solution in pH 3.56 bitartrate buffer was sampled into the syringe and maintained at the room temperature for 1, 5, or 10 min. Replicate samples (n = 5) were run for each equilibration time using 50 mL headspace samples.



**Determination of Henry's coefficient of SO<sub>2</sub> as a function of ethanol concentration.** SO<sub>2</sub> solutions (1 mg/L) were prepared with varying ethanol concentrations (0, 8, 11, 14, or 17% v/v) at pH 0 by adding SO<sub>2</sub> stock directly into a syringe containing 10 mL of appropriate ethanol/HCl solutions. At this pH, SO<sub>2</sub> was predominantly in the molecular SO<sub>2</sub> form and care was taken to keep the syringe closed except during sampling. Replicate samples (n = 8) were analyzed by the HS-GDT method for each ethanol concentration. The Henry's coefficient (K<sub>H</sub>) of SO<sub>2</sub> was then calculated in units of Atm L/mol for each ethanol concentration as  $K_H = P_{SO_2} / [SO_2]$ .

**Determination of pK<sub>a</sub> of SO<sub>2</sub> as a function of ethanol.** SO<sub>2</sub> solutions (50 mg/L) were prepared in a potassium bitartrate buffer with one of four ethanol concentrations (0, 7, 14, or 20% v/v) at one of four pH values (3.1, 3.3, 3.5, or 3.7), for a total of 16 solutions. pH was adjusted by dropwise addition of either 25% phosphoric acid or 1 N sodium hydroxide. The HS-GDT protocol was then used to determine P<sub>SO<sub>2</sub></sub>, which was then converted to molecular SO<sub>2</sub>. The apparent pK<sub>a</sub> for each solution was calculated from the measured molecular SO<sub>2</sub> and known free SO<sub>2</sub> concentration using the Henderson-Hasselbalch equation. The room temperature at the time of the experiment was 21.9°C.

**Figures of merit for the HS-GDT method: linearity, % CV, and detection limit.** Standard solutions at nominal concentrations of 5, 10, 15, 20, 30, and 40 mg/L SO<sub>2</sub> were made in potassium bitartrate buffer (pH 3.56) with 12% (v/v) ethanol and at nominal concentrations of 5, 10, 25, 50, and 75 mg/L SO<sub>2</sub> in aqueous potassium bitartrate buffer (pH 3.56). Actual free SO<sub>2</sub> concentrations were determined by the Ripper method (Iland et al. 2004). Each sample was analyzed by the HS-GDT method in replicate (n = 12; 72 analyses total). The coefficient of variation for each standard concentration was calculated as the standard deviation divided by the mean and expressed as a percentage. Linearity was evaluated by linear regression of measured P<sub>SO<sub>2</sub></sub> versus expected P<sub>SO<sub>2</sub></sub>. The detection limit (LOD) was calculated as  $3 \times s$ , where s is the signal-independent noise calculated by Pallesen's method (Berthouex and Brown 2002).

**Comparison of A-O and HS-GDT for measuring SO<sub>2</sub> in commercial wines.** Twenty-seven commercial wines (9 red, 14 white, and 4 blush) were evaluated by both the A-O and HS-GDT methods on site. The red wines were from the 2009 to 2012 vintages and the white wines were from the 2011 to 2013 vintages. The wines were from Australia, Argentina, and the United States (California, New York, Oregon, and Washington). Composition of the commercial wines was determined by ETS Laboratories (St. Helena, CA) using accredited methods for pH (by meter), alcohol (FTIR), and glucose plus fructose (enzymatic); additionally tannins, polymeric, and total anthocyanins were measured by HPLC (Waterhouse et al. 1999). Similar to standard curves, the LOD for molecular SO<sub>2</sub> was calculated as  $3 \times s$ , where s is the signal-independent noise calculated by Pallesen's method (Berthouex and Brown 2002).

To evaluate whether discrepancies between A-O and HS-GDT measurements could be explained by the presence of

labile anthocyanin-bisulfite adducts, the concentration of anthocyanin-bisulfite adducts, [Flav-Bisulfite], was estimated by solving a series of equilibrium equations (Equations 5 to 8):

$$K_D = [\text{Flav}] * [\text{Bisulfite}_{\text{HS}}] / [\text{Flav-Bisulfite}] \quad \text{Eq. 5}$$

$$K_A = [\text{Quin}][\text{H}^+] / [\text{Flav}] \quad \text{Eq. 6}$$

$$[\text{Total Anth}] = [\text{Flav}] + [\text{Quin}] + [\text{Flav-Bisulfite}] \quad \text{Eq. 7}$$

$$[\text{Bisulfite}_{\text{AO}}] = [\text{Bisulfite}_{\text{HS}}] + [\text{Flav-Bisulfite}] \quad \text{Eq. 8}$$

[Flav] and [Quin] are the concentrations of flavylum and quinoid forms of anthocyanins. [Flav-Bisulfite] is the concentration of anthocyanin bisulfite adducts, and the free SO<sub>2</sub> concentration measured by A-O, [Bisulfite<sub>AO</sub>], was assumed to equal the sum of [Flav-Bisulfite] and [Bisulfite<sub>HS</sub>], i.e., all of the difference between A-O and HS-GST could be explained by binding to flavylum ions. Acid dissociation constants (pK<sub>a</sub>) for anthocyanins were based on literature values and assumed to be 2.94 for *Vitis vinifera* wines and 2.6 for diglucoside-containing hybrid-based wines (Brouillard and El Hage Chahine 1980, Francis and Markakis 1989, Ribéreau-Gayon et al. 2007). The dissociation constant (K<sub>d</sub>) for Flav-Bisulfite adducts was assumed to be  $1 \times 10^{-5}$  (Burrroughs 1975, Timberlake and Bridle 1976). [Total Anth] was the total anthocyanin concentration, and [H<sup>+</sup>] for each wine was derived from the measured pH. For each wine, K<sub>D</sub>, K<sub>A</sub>, [Total Anth], [H<sup>+</sup>], and [Bisulfite<sub>AO</sub>] were independent variables and the equations were solved for [Flav], [Quin], [Flav-Bisulfite], and [Bisulfite<sub>HS</sub>] using R v.3.1.2 (The R Foundation for Statistical Computing, Vienna, Austria) and the "nleqslv" v.2.5 add-in package.

**Statistical Analyses.** JMP Pro 11.0.0 (SAS Institute, Inc., Cary, NC) was used for statistical analyses. Comparison of means tests for the effects of ethanol on pK<sub>a</sub> and K<sub>H</sub> were performed by one-way ANOVA followed by a Tukey test if the effect was significant. Linear regression analyses for calibration curve comparison of A-O versus HS-GDT data were performed with 1/x weighting factors. Linear regressions of wine composition parameters and the difference between A-O and HS-GDT SO<sub>2</sub> values were performed without weighting.

## Results and Discussion

**Apparatus and materials for SO<sub>2</sub> measurements by HS-GDT.** The apparatus was constructed and used in the protocol as described (Figure 1). An example of a partially discolored Gastec 5Lb tube is also shown (Supplemental Figure 1). Three different gas detector tubes were initially investigated: i) Gastec 5Lb, which rely on the reaction of SO<sub>2</sub> with BaCl<sub>2</sub> to generate HCl, resulting in color change of a pH-sensitive dye; ii) Sensidyne/Kitigawa 103SE SO<sub>2</sub> tubes, which rely on the colorimetric reaction of SO<sub>2</sub> with NaOH to generate Na<sub>2</sub>SO<sub>3</sub>; and iii) Draeger 0.1/a tubes, which rely on reaction of SO<sub>2</sub> with NaHgCl<sub>4</sub> to generate HCl, whose evolution is detected by a pH-sensitive indicator. Initial trials with alcoholic buffer solutions containing SO<sub>2</sub> resulted in no detectable signal for the Draeger 0.1/a tubes, and thus these were not further considered. We observed inconsistent results with Sensidyne

tubes, with coefficients of variance (CV) >20% in the length of tube darkened, while Gastec 5Lb tubes yielded reproducible results (CV <10%). The reason for the poor performance of some of the tubes is unclear. However, commercial tubes require that gas sampling occur at precise, non-linear, and proprietary rates that are unique to each manufacturer (Haag 2001) and the poor performance may reflect incompatibility of some brands of tubes with this method's syringe-depression sampling rate. The Gastec 5Lb tubes were selected for further study, as they performed best with the current syringe apparatus sampling protocol.

**Evaluation of interferences.** Gastec GDT exhibited uniform discoloration from green to light blue even when exposed to “blank” 14% ethanol buffer solutions (data not shown). However, because this color change was distinct from the color change associated with SO<sub>2</sub> (green to yellow), it did not interfere with measurements. Because the reaction chemistry of the GDT employed will respond to all volatile acids, we investigated the susceptibility of the HS-GDT method to interference from acetic (1.4 g/L), lactic (5 g/L), and sorbic (300 mg/L) acids. A small signal, equivalent to <0.1 mg/L molecular SO<sub>2</sub>, was observed at these concentrations (data not shown). Since these concentrations represent the upper end of concentrations expected to be found in wines (Zoecklein et al. 1995), the susceptibility of the method to these volatile acid interferences is low.


**Interpolation of manufacturer markings on GDT.** The manufacturer-provided markings on the GDT were used to quantify the SO<sub>2</sub> vapor pressure, P<sub>SO<sub>2</sub></sub>. To facilitate interpolation of markings, best-fit cubic functions relating P<sub>SO<sub>2</sub></sub> to distance along the tube were generated for each lot of Gastec GDT (Supplemental Figure 2, Supplemental Table 1). P<sub>SO<sub>2</sub></sub> for each sample could then be calculated using Equation 1, which in turn could be used to calculate molecular SO<sub>2</sub> and free SO<sub>2</sub>. Tube-to-tube variation in marker spacing within the same lot number was imperceptible and not a source of error.

**Effect of equilibrium time on HS-GDT measurements.** Fifty mg/L SO<sub>2</sub> solutions in pH 3.56 buffer were analyzed by the HS-GDT apparatus using varying static equilibration times (1, 5, or 10 min) at room temperature. No significant differences were observed, indicating that headspace and liquid-phase SO<sub>2</sub> equilibrate quickly. Five minutes was selected as the optimal equilibrium time for the HS-GDT method because it yielded the best precision (Table 1).

**Effects of ethanol on Henry's coefficient of SO<sub>2</sub>.** Using the HS-GDT assay, Henry's coefficient (K<sub>H</sub>) was determined for 1 mg/L SO<sub>2</sub> solutions at 21°C over varying ethanol con-

centrations. Samples were acidified with HCl to ~pH 0 to favor the molecular SO<sub>2</sub> form and eliminate any confounding ethanol-induced effects on the pK<sub>a</sub> of SO<sub>2</sub>. K<sub>H</sub> was independent of ethanol concentration over the range of 0 to 17% v/v ethanol (Figure 2) and was calculated to be 0.38 Atm/M at 21°C, comparable to a recent previous report of 0.28 Atm/M for SO<sub>2</sub> in water at 25°C (Zhang et al. 2013). The lack of dependence of SO<sub>2</sub> volatility on ethanol concentration is in contrast to other wine volatiles. For example, the volatility of many wine esters is about two-fold lower in model wine than in a model juice solution (Mouret et al. 2012). The minimal effect of ethanol on SO<sub>2</sub> volatility may result from its greater polarity than many other wine volatiles. Regardless of the explanation, these results indicate that it is not necessary to control for ethanol concentration prior to HS-GDT analyses of wine.

**Effect of ethanol concentration on the pK<sub>a</sub> of SO<sub>2</sub>.** P<sub>SO<sub>2</sub></sub> was determined by HS-GDT for model wine solutions with different ethanol concentrations and pH values and the P<sub>SO<sub>2</sub></sub> and pH values were subsequently used to calculate pK<sub>a</sub> as a function of ethanol concentration (Figure 3). Because K<sub>H</sub> is independent of ethanol concentration, changes in P<sub>SO<sub>2</sub></sub> could be assigned solely to effects on pK<sub>a</sub>. The best fit line ( $y = 0.0137x + 1.83$ ,  $r^2 = 0.97$ ) yielded predicted acid-dissociation constants of pK<sub>a</sub> = 1.83 at 0% ethanol and pK<sub>a</sub> = 1.96 at 10% alcohol. These values are comparable to pK<sub>a</sub> values determined by titrimetry in model solutions for 0 and 10% ethanol solutions (1.81 and 2.00, respectively; Table 2; Usseglio-Tomasset and Bosia 1984). Interestingly, the pK<sub>a</sub> of SO<sub>2</sub> in water (1.81) is widely recommended in calculations of molecular SO<sub>2</sub> from free SO<sub>2</sub> in wine regardless of the ethanol concentration, which is expected to lead to 25 to 50% underestimates of the actual molecular SO<sub>2</sub> concentration in wines, depending on the pH. Temperature and ionic strength can also impact pK<sub>a</sub> to a lesser degree, and should be taken into consideration (Usseglio-Tomasset and Bosia 1984). In this paper, all conversions between free and molecular SO<sub>2</sub> are based on ethanol-corrected pK<sub>a</sub> values at specific room temperatures and 50 mM ionic strength unless otherwise specified.




**Figure 2** SO<sub>2</sub> vapor pressure as a function of ethanol concentration. Samples were prepared at ~pH 0. Error bars represent standard errors (n = 8).

**Table 1** Effect of equilibrium time on SO<sub>2</sub> vapor pressure measured by headspace gas detection tube.

	Equilibration time		
	1 min	5 min	10 min
P <sub>SO<sub>2</sub></sub> (µL/L)	1.18	1.19	1.21
Standard deviation	0.16	0.11	0.23
Coefficient of variation (%)	13.28	9.37	19.29

**Figures of merit for the HS-GDT method.** Calibration curves were analyzed by the HS-GDT method using SO<sub>2</sub> standards prepared at pH 3.56 in either aqueous buffer (0.09 to 1.28 mg/L molecular SO<sub>2</sub>) or model wine (0.14 to 1.12 mg/L molecular SO<sub>2</sub>) using either a 50 mL or 200 mL (4 × 50 mL) headspace sample. Limits of detection and %CV for commercial wines were also calculated (Table 3). The slopes and intercepts for all methods were nearly identical (<5% difference for each parameter across the four curves), indicating



**Figure 3** Experimentally determined  $pK_1$  values for SO<sub>2</sub> in a buffered standard solution as a function of alcohol concentration. Error bars represent standard deviations ( $n = 4$ ).

**Table 2** Comparison of experimentally determined  $pK_1$  values to literature  $pK_1$  values for SO<sub>2</sub>.

Alcohol	$pK_1$ (observed) <sup>a</sup>	$pK_1$ (literature) <sup>b</sup>
0%	1.83	1.80
5%	1.90	1.90
10%	1.96	2.00
15%	2.03	2.10
20%	2.10	2.20

<sup>a</sup>Experimentally derived,  $T = 21.8^\circ\text{C}$ .

<sup>b</sup>Calculated based on formulae from Usseglio-Tomasset and Bosia (1984) using  $T = 21.8^\circ\text{C}$  and ionic strength = 0.056M.

that the method was scalable with headspace volume size and that using Equation 1 was valid. As expected from our observations that wine-like concentrations of ethanol do not affect Henry's coefficient, the best fit lines for the aqueous and model wine calibrations were nearly identical ( $y = 5.92x - 0.48$  versus  $y = 5.99x - 0.55$  for 200 mL sample sizes).

Although reasonable detection limits (0.21 mg/L molecular SO<sub>2</sub>) could be achieved for 50 mL sample sizes in the aqueous standards, detection limits were considerably higher for model wine (0.59 mg/L molecular SO<sub>2</sub>). Additionally, linearity was unacceptable for the 50 mL samples from model wine ( $r^2 = .75$ ), and precision over the range of molecular SO<sub>2</sub> typically targeted in wine (>0.4 mg/L) was also mediocre (mean CV = 18%). In the model wine system, using repeat sampling improved limits of detection from 0.59 mg/L for a 50 mL headspace sample to 0.21 mg/L for a 200 mL sample. Repeat sampling also improved the precision (mean CV = 8% for standards with >0.4 mg/L molecular SO<sub>2</sub>). Based on these results, 200 mL headspace sample sizes (i.e., repeated sampling) were used for later analyses on commercial wines.

Unexpectedly, the y-intercept was negative for all calibration curves, indicating that a minimal amount of SO<sub>2</sub> (~0.07 mg/L molecular SO<sub>2</sub>) was necessary to effect a color change. An explanation for this phenomenon is not clear, but may be an effect of the differences between the manufacturer's optimized flow rate through the GDT and the flow rate used in our experiments (Haag 2001). Regardless, the offset is well below typical recommendations for molecular SO<sub>2</sub> necessary to prevent spoilage, and thus the method should still be useful for most wines.

**Comparison of HS-GDT and A-O methods.** Twenty-seven commercial wines (9 red, 14 white, and 4 blush) were evaluated by both HS-GDT and A-O methods. HS-GDT measurements in  $P_{\text{SO}_2}$  were converted to free SO<sub>2</sub> values. The wines had a wide range of free SO<sub>2</sub> concentrations (4 to 51 mg/L by A-O). The average precision for HS-GDT measurements of wines with >0.4 mg/L molecular SO<sub>2</sub> (mean CV = 17%) was different than values achieved with calibration standards, and the limit of detection (0.06 mg/L) was three-fold lower than the values achieved for calibration standards (Table 3). The reason for the reduced precision in real wines may be due to limiting the number of repetitions to those that

**Table 3** Figures of merit for the headspace gas detection tube (HS-GDT) method. Aqueous and 10% ethanol model wines both at pH 3.56; commercial wines are described in the text.

HS-GDT conditions	Intercept <sup>a</sup>	Slope	$r^2$	Mean CV for molecular SO <sub>2</sub> >0.4 mg/L (%) <sup>b</sup>	Detection limit for molecular SO <sub>2</sub> (mg/L) <sup>c</sup>
50 mL headspace sample, aqueous <sup>d</sup>	-0.48	5.84	0.99	15%	0.21
50 mL headspace sample, model wine	-0.60	5.96	0.75	18%	0.59
200 mL headspace sample, aqueous	-0.48	5.92	0.97	9%	0.29
200 mL headspace sample, model wine	-0.55	5.99	0.99	8%	0.21
200 mL headspace sample, commercial wines	NA <sup>e</sup>	NA	NA	17%	0.06

<sup>a</sup>Intercept, slope, and  $r^2$  are based on  $1/x$  weighted best-fit lines for plots of  $P_{\text{SO}_2}$  ( $\mu\text{L/L}$ ) versus molecular SO<sub>2</sub> (mg/L).

<sup>b</sup>CV was calculated as standard deviation/mean.

<sup>c</sup>Detection limit was calculated as  $3 \times$  signal-independent noise,  $\sigma$ .

<sup>d</sup>Responses for 50 mL sample sizes were scaled by a factor of four, as described in Materials and Methods.

<sup>e</sup>NA: Not applicable.




would fit on one tube, thus increasing the variability on the wines with higher concentrations of SO<sub>2</sub>. The lower limit of detection for HS-GDT in real wines as opposed to calibration standards is also unclear. Potentially, the imprecision observed in model wine for low concentration standards (<0.2 mg/L) arises in part from variable degrees of oxidation. In real wines, these losses may be buffered by dissociation of the bound SO<sub>2</sub> pool or minimized by the presence of the other antioxidant wine components.

A very good correlation was observed between the HS-GDT and A-O methods in white and blush wines ( $r^2 = 0.97$ ), but much weaker correlation was observed for red wines ( $r^2 = 0.72$ , Figure 4). HS-GDT values for red wines were on average only 49% (range = 24 to 76%) of the A-O value recorded for the same wine, with absolute differences approaching 25 mg/L for some wines. This large discrepancy between the two methods for red wines is likely an artifact of A-O sample preparation steps (acidification, dilution) and long sampling time (10 to 15 min) that allow dissociation of weakly-bound anthocyanin-bisulfite adducts (Burroughs 1975) and erroneously high measurements of free SO<sub>2</sub>. Although equilibrating the 10 mL sample by contact with 50 mL headspace in the HS-GDT method will slightly decrease wine SO<sub>2</sub> concentrations due to volatilization, the effect will be minimal (<1% change in bisulfite or molecular SO<sub>2</sub> at equilibrium) because the pH does not change. Similar results have been observed in the limited number of other studies that have used non-perturbing methods to quantify free SO<sub>2</sub> in red wines. For example, free SO<sub>2</sub> values were up to an order of magnitude higher by Ripper and A-O than when measured by CE, where equilibria were not perturbed (Bogren 1996). A group using a non-perturbing headspace GC method reported a 45% higher free SO<sub>2</sub> value in a red wine than with A-O, although this

report used only a single wine adjusted to different pH values rather than multiple wines (Davis et al. 1983). The HS-GDT and A-O methods were in better agreement for white and blush wines. Based on the slope of the regression analysis (Figure 4), free SO<sub>2</sub> values by HS-GDT were 87% that of the A-O method. This slight discrepancy (13% lower by HS-GDT) is comparable to that observed in a previous report comparing CE and A-O (Bogren 1996). This previous work suggested that any discrepancy could be explained by dissociation of weakly bound adducts of bisulfite, diacetyl, pyruvate, and other carbonyl species during A-O analyses. A separate report comparing HS-GC and A-O methods for a single white wine reported differences of <5% (Davis et al. 1983), but this work used a pK<sub>a</sub> value for SO<sub>2</sub> (1.81) that was uncorrected for the effects of ethanol, which would have resulted in a higher free SO<sub>2</sub> value based on measured headspace SO<sub>2</sub>. In contrast, a good correlation was observed between A-O and a modified A-O approach using GDT for both red and white wines (Pegram et al. 2013); in that approach, the wine pH was decreased to <2 to favor SO<sub>2</sub> volatilization, resulting in perturbation of equilibrium conditions and measurement of weak bisulfite adducts.

**Predicted versus actual discrepancies in molecular SO<sub>2</sub> in red wine: the role of anthocyanin bisulfite adducts.** To evaluate the hypothesis that the large differences between A-O and HS-GDT methods for red wines could be explained by dissociation of anthocyanin-bisulfite adducts, the concentration of anthocyanin-bisulfite adducts was estimated by solving a system of nonlinear equations involving the equilibria of true free HSO<sub>3</sub><sup>-</sup> and flavylium, quinoid, and bisulfite adduct forms of anthocyanins. Measured values of total anthocyanins, pH, and free SO<sub>2</sub> by A-O (assumed to be true free SO<sub>2</sub> + anthocyanin-bound adducts) were used in these calculations; a reported literature value of 2.94 for the anthocyanin hydration constant (pK<sub>h</sub>) was used for anthocyanins from *V. vinifera* wines (Brouillard and El Hage Chahine 1980) and a value of pK<sub>h</sub> = 2.6 was used for the diglycoside-containing, hybrid-based wines (Francis and Markakis 1989). Only adducts of flavylium ion forms and bisulfite were considered in the model, as the other anthocyanin forms do not bind bisulfite (Timberlake and Bridle 1967). The effects of self-aggregation on flavylium cations were ignored in the model. The apparent equilibrium constant for dimerization, K<sub>d</sub>, was calculated to be no greater than 6000 for all red wine samples (estimated from Eq A11 in Houbiers et al. 1998), which would have resulted in only a small overestimation (less than 10%) of the flavylium ion concentrations used in our model. The observed discrepancy in molecular SO<sub>2</sub> between the A-O and HS-GDT methods was plotted against the estimated concentration of anthocyanin-bisulfite adducts, expressed as SO<sub>2</sub> equivalents (Figure 5A); it exhibits a strong positive correlation ( $r^2 = 0.936$ ) and has a slope near unity (0.94), supporting our hypothesis. Similarly, molecular SO<sub>2</sub> by HS-GDT plotted against A-O values that were corrected for the calculated dissociation (Figure 5B) yielded a modest  $r^2$  of 0.755 and a slope of 0.765, very similar to the curve for the white and blush wines in Figure 4. Based on linear regression



**Figure 4** Plot of molecular SO<sub>2</sub> measured by the headspace gas detection tube (HS-GDT) method versus an aeration-oxidation (A-O) reference method for 9 red, 14 white, and 4 blush wines. Each data point represents the SO<sub>2</sub> concentration of an individual wine (mean of n = 3). Best-fit line equations determined with 1/x weighting were:  $y = 0.82x + 0.05$ ,  $r^2 = 0.97$  (white and blush wines, open circles) and  $y = 0.32x + 0.05$ ,  $r^2 = 0.72$  (red wines, solid circles).

analyses, other measured wine components (tannins, glucose + fructose, ethanol) did not correlate with the percentage difference between A-O and HS-GDT methods in either red or white wines (data not shown). The lack of correlation of method difference with glucose concentration is not surprising. Dissociation of glucose-bisulfite complexes has been suggested as a potential problem to standard approaches to SO<sub>2</sub> measurement (Burroughs and Sparks 1964) and several wines included in the study had fructose + glucose concentrations >50 g/L. However, the first order rate constant for dissociation of the glucose-bisulfite complex is reportedly  $3.7 \times 10^{-4}/\text{min}$  at ~pH 1 (Vas 1949), which is slow compared to the time necessary for A-O or Ripper analysis and considerably slower than the first order rate constant for dissociation of anthocyanin-bisulfite adducts (0.2/min) (Brouillard and El Hage Chahine 1980).


## Conclusion

The HS-GDT method described here provides a convenient and inexpensive means to measure molecular and free SO<sub>2</sub> concentrations typically found in wines without perturbing the equilibria of free and bound SO<sub>2</sub> forms. Our work also supports previous (and often overlooked) observations that standard approaches to SO<sub>2</sub> measurement overestimate free and molecular SO<sub>2</sub> due to dissociation of weakly bound bisulfite adducts. This is particularly notable in red wines, where discrepancies up to five-fold were noted, and is also a minor issue in white and blush wines. Thus, the HS-GDT method may be of use to researchers interested in determining the mechanisms of wine oxidation or better establishing microbial tolerance to SO<sub>2</sub>. Because of the low requirements for consumables and fast analysis time (~5 min), HS-GDT could also be adopted for use by winemakers to quantify SO<sub>2</sub> without a specialized lab space.

Although additional validation of the HS-GDT method is desirable, traditional approaches to method validation, such as standard addition and percent recovery calculations, are particularly challenging when working with free (or molecular) SO<sub>2</sub> in wine matrices. Unlike the validation of a total SO<sub>2</sub> method, when any addition of SO<sub>2</sub> to a wine matrix would be reflected in the same increase in total SO<sub>2</sub>, addition of SO<sub>2</sub> will result in variable increases in free or molecular SO<sub>2</sub> because of varying concentrations of SO<sub>2</sub> binders across wines. Standard addition could potentially be done by measuring major SO<sub>2</sub> binders in wine to calculate the expected increase in free SO<sub>2</sub> following SO<sub>2</sub> addition. Because of the challenges of using standard addition, further validation of the HS-GDT method could be done by comparison with free (or molecular) SO<sub>2</sub> values determined by other non-perturbing techniques, such as headspace GC, headspace IR, and CE.

## Literature Cited

- AOAC. 2007. Official Methods of Analysis of AOAC International, 18th ed., Horwitz (ed.). Online. AOAC International, Gaithersburg, MD.
- AOAC. 2015. Official Methods of Analysis of AOAC International, 20th ed., Horwitz (ed.). Online. AOAC International, Gaithersburg, MD.
- Beech, F.W., L.F. Burroughs, C.F. Timberlake, and G.C. Whiting. 1979. Recent progress on chemical aspects and the antimicrobial action of SO<sub>2</sub>. *Bull. OIV* 52:1001-1022.
- Berthouex, P.M., and L.C. Brown. 2002. *Statistics for Environmental Engineers*. Lewis Publishers, Boca Raton, FL.
- Bogren, L.N.H. 1996. The effect of sulfur dioxide on *Saccharomyces cerevisiae*. Thesis, UC Davis, Davis, CA.
- Boulton, R.B., V.L. Singleton, L.F. Bisson, and R.E. Kunkee. 1996. *Principles and Practices of Winemaking*. Chapman & Hall, New York.
- Brouillard, R., and J.M. El Hage Chahine. 1980. Chemistry of anthocyanin pigments. 6. Kinetic and thermodynamic study of hydrogen sulfite addition to cyanin. Formation of a highly stable Meisenheimer-type adduct derived from a 2-phenylbenzopyrylium salt. *J. Am. Chem. Soc.* 102:5375-5378.
- Buechsenstein, J.W., and C.S. Ough. 1978. SO<sub>2</sub> determination by aeration-oxidation: A comparison with Ripper. *Am. J. Enol. Vitic.* 29:161-164.
- Burroughs, L.F. 1975. Determining free sulfur dioxide in red wine. *Am. J. Enol. Vitic.* 26:25-29.



**Figure 5** (A) Plot of differences in molecular SO<sub>2</sub> by aeration-oxidation (A-O) and headspace gas detection tube (HS-GDT) methods versus estimated decrease in molecular SO<sub>2</sub> resulting from anthocyanin-bisulfite adduct formation. Best fit-lines equations determined with 1/x weighting were  $y = 0.9436x - 0.0161$ ,  $r^2 = 0.9362$ . (B) Plot of molecular SO<sub>2</sub> measured by the HS-GDT method versus the A-O reference method, corrected for theoretical anthocyanin-bisulfite adduct formation for nine red wines. Each data point represents the SO<sub>2</sub> concentration of an individual wine (mean of n = 3). Best fit-lines equation determined with 1/x weighting were:  $y = 0.7652x + 0.0461$ ,  $r^2 = 0.76$ .



- Burroughs, L.F., and A.H. Sparks. 1964. The determination of the free sulphur dioxide content of ciders. *Analyst* 89:55-60.
- Burroughs, L.F., and A.H. Sparks. 1973. Sulphite binding power of wines and ciders. I. equilibrium constants for dissociation of carbonyl bisulphite compounds. *J. Sci. Food Agric.* 24:187-198.
- Čmelík, J., J. Machát, E. Niedobová, V. Otruba, and V. Kanický. 2005. Determination of free and total sulfur dioxide in wine samples by vapour-generation inductively coupled plasma-optical-emission spectrometry. *Anal. Bioanal. Chem.* 383:483-488.
- Collins, T., and R.B. Boulton. 1996. The analysis of free sulfur dioxide and sugars in juices and wines by capillary electrophoresis. *In* *Proceeding of Oenologie 95: 5th Symposium International d'Oenologie*. A. Lonvaud-Funel (ed.), pp. 637-640. Lavoisier, Paris, France.
- Danilewicz, J.C. 2011. Mechanism of autoxidation of polyphenols and participation of sulfite in wine: Key role of iron. *Am. J. Enol. Vitic.* 62:319-328.
- Davis, E.G., D. Barnett, and P.M. Moy. 1983. Determination of molecular and free sulfur dioxide in foods by headspace gas chromatography. *J. Food Technol.* 18:233-240.
- Divol, B., M. du Toit, and E. Duckitt. 2012. Surviving in the presence of sulphur dioxide: Strategies developed by wine yeasts. *Appl. Microbiol. Biol.* 95:601-613.
- Francis, F.J., and P.C. Markakis. 1989. Food colorants: Anthocyanins. *Crit. Rev. Food Sci.* 28:273-314.
- Godden, P., L. Francis, J. Field, M. Gishen, A. Coulter, P. Valente, P. Høj, and E. Robinson. 2001. Wine bottle closures: Physical characteristics and effect on composition and sensory properties of a Semillon wine. I. Performance up to 20 months post-bottling. *Aust. J. Grape Wine Res.* 7:64-105.
- Haag, W.R. 2001. Interchangeability of gas detection tubes and hand pumps. *Am. Ind. Hyg. Assoc. J.* 62:65-69.
- Henningsen, J., and J. Hald. 2003. Quantitative analysis of dilute mixtures of SO<sub>2</sub> in N<sub>2</sub> at 7.4 μm by difference frequency spectroscopy. *Appl. Phys. B-Lasers O.* 76:441-449.
- Houbiers, C., J.C. Lima, A.L. Maçanita, and H. Santos. 1998. Color stabilization of malvidin 3-glucoside: Self aggregation of the flavylum cation and copigmentation with the Z-chalcone form. *J. Phys. Chem. B.* 102:3578-3585.
- Huang, M.D., H. Becker-Ross, S. Florek, U. Heitmann, M. Okrus, and C.D. Patz. 2008. Determination of sulfur forms in wine including free and total sulfur dioxide based on molecular absorption of carbon monosulfide in the air-acetylene flame. *Anal. Bioanal. Chem.* 390:361-367.
- Iland, P., N. Bruer, G. Edwards, S. Weeks, and E. Wilkes. 2004. *Chemical Analysis of Grapes and Wine: Techniques and Concepts*. Patrick Iland Wine Promotions, Campbelltown, S. Australia.
- Joslyn, M.A. 1955. Determination of free and total sulfur dioxide in white table wines. *J. Agric. Food Chem.* 3:686-695.
- Mouret, J.R., S. Morakul, P. Nicolle, V. Athes, and J.M. Sablayrolles. 2012. Gas-liquid transfer of aroma compounds during winemaking fermentations. *In* *LWT-Food Science and Technology* 49: Euro-Mediterranean Symposium on Fruit and Vegetable Processing. C. Renard (ed.), pp. 238-244. Elsevier Ltd, Waltham, MA.
- Park, S.K. 2008. Development of a method to measure hydrogen sulfide in wine fermentation. *J. Microbiol. Biotechnol.* 18:1550-1554.
- Pegram, Z., M.T. Kwasniewski, and G.L. Sacks. 2013. Simplified method for free SO<sub>2</sub> measurement using gas detection tubes. *Am. J. Enol. Vitic.* 64:405-410.
- Pundir, C.S., and R. Rawal. 2013. Determination of sulfite with emphasis on biosensing methods: A review. *Anal. Bioanal. Chem.* 405:3049-3062.
- Rankine, B.C., and K.F. Pocock. 1970. Alkalimetric determination of sulphur dioxide in wine. *Aust. Wine Brew. Spirit Rev.* 88:40-44.
- Ribéreau-Gayon, P., D. Dubourdieu, B. Donèche, and A. Lonvaud. 2007. *Handbook of Enology*. Vol. 1. John Wiley, Hoboken, NJ.
- Ruiz-Capillas, C., and F. Jiménez-Colmenero. 2009. Application of flow injection analysis for determining sulphites in food and beverages: A review. *Food Chem.* 112:487-493.
- Sander, R. 2015. Henry's law constants. *In* *NIST Chemistry WebBook, NIST Standard Reference Database Number 69*, P.J. Linstrom and W.G. Mallard (eds.). National Institute of Standards and Technology, Gaithersburg MD. <http://webbook.nist.gov> (accessed 14 Feb 2015).
- Timberlake, C.F., and P. Bridle. 1967. Flavylum salts, anthocyanidins and anthocyanins. II.—Reactions with sulphur dioxide. *J. Sci. Food Agric.* 18:479-485.
- Timberlake, C.F., and P. Bridle. 1976. Interactions between anthocyanins, phenolic compounds, and acetaldehyde and their significance in red wines. *Am. J. Enol. Vitic.* 27:97-105.
- Ugliano, M., and P.A. Henschke. 2010. Comparison of three methods for accurate quantification of hydrogen sulfide during fermentation. *Anal. Chim. Acta* 660:87-91.
- Usseglio-Tomasset, L., and P.D. Bosia. 1984. La prima costante di dissociazione dell'acido solforoso. *Vini d'Italia* 26:7-14.
- Vas, K. 1949. The equilibrium between glucose and sulphurous acid. *J. Chem. Technol. Biotechnol.* 68:340-343.
- Waterhouse, A.L., and R.J. Elias. 2010. Chemical and physical deterioration of wine. *In* *Chemical Deterioration and Physical Instability of Food and Beverages*. L.H. Skibsted et al. (eds.), pp. 446-482. Woodhead Publishing, Boca Raton, FL.
- Waterhouse, A.L., S.F. Price, and J.D. McCord. 1999. Reversed-phase high-performance liquid chromatography methods for analysis of wine polyphenols. *In* *Methods of Enzymology*, Vol. 299: Oxidants and Antioxidants Part A. L. Packer (ed.), pp. 113-121. Elsevier, Waltham, MA.
- Zhang, N., J.B. Zhang, Y.F. Zhang, J. Bai, and X.H. Wei. 2013. Solubility and Henry's law constant of sulfur dioxide in aqueous polyethylene glycol 300 solution at different temperatures and pressures. *Fluid Phase Equilib.* 348:9-16.
- Zoecklein, B.W., K.C. Fugelsang, B.H. Gump, and F.S. Nury. 1995. *Wine Analysis and Production*. Springer, New York.